



Mortuary - Autopsy Room

Room Code	MOR-AU
Briefed Area	30.00 m ²
Ceiling Height	3.0 m
Occupancy	up to 4 staff
Hours of Operation	8am to 5pm, daily (may vary depending on services requirements)

Description | The Mortuary - Autopsy Room is for undertaking post mortem examinations. An observation area is to be provided to at least one Autopsy Room within a department.

Electrical | PROTECTION: body protected ☐
PROTECTION: cardiac protected ☐

Lighting | LIGHTING: general ☐
LIGHTING: colour corrected ☒
LIGHTING: dimmable ☐
LIGHTING: indirect ☐

Nurse Call and Duress | NURSE CALL SYSTEM: buttons / handset ☐

HVAC | AIRCONDITIONING: general ☐
AIRCONDITIONING: HEPA filtered ☐
AIRCONDITIONING: positive pressure ☐
AIRCONDITIONING: negative pressure ☒
VENTILATION: exhaust ☒
VENTILATION: natural ☐

Additional Considerations |

- Layout of tables and orientation of observation bay may be adjusted to suit the type of autopsy being performed (e.g. coronial or non-coronial). The Mortuary - Autopsy Room should be located with ready access to Mortuary - Body Holding.
- doors to be sealed for odours, post mortem table to have water flushing system. Appropriate Smoke/Heat detection to be considered.



Doors and Windows

CODE	DESCRIPTION	COMMENT
DOHI-114.01	DOOR: hinged, 1 1/2 leaves, 1400 clear opening, solid, standard vision panel	[DWGL-033] to corridor; door seals required
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to corridor side, extent dependent on movement of beds and mobile equipment in adjacent area; to be coordinated with corridor wall protection
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to room side
AFDPR-056.03	DOOR FRAME PROTECTION: full wrap, to 1200H	[DWPR-010] optional, provision and extent dependent on frame material/finish and movement of beds and mobile equipment in adjacent area
DOHI-013.01	DOOR: hinged, 1 leaf, 900 clear opening, solid, standard vision panel	[DWGL-005] to anteroom; door seals required
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to room side
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to corridor side, extent dependent on movement of beds and mobile equipment in adjacent area; to be coordinated with corridor wall protection
WIFX-008.01	WINDOW: fixed, internal, single glazed, sill at 1050H	[DWOB-005] to observation area if provided



Finishes

CODE	DESCRIPTION	COMMENT
FLVY-112	FLOOR FINISH: vinyl, seamless, non-slip, safety	[FLVY-004] junction with evisceration table to be considered to prevent water ingress
FLSK-021	SKIRTING: vinyl, integral with floor vinyl, coved	[FLSK-010]
WLFI-001	WALL FINISH: paint	[WLWA-004]
WLFI-011.06	WALL FINISH: vinyl, to 2100 AFFL	[WLWA-023]
CLFS-011	CEILING: flush set, suspended	[CLFS-006]
CLFI-002	CEILING FINISH: paint, clinical areas	[CLFS-006]
CLCN-031	CORNICE: square set	[CLCN-008]



Joinery

CODE	DESCRIPTION	GROUP	QTY	COMMENT
JOBE-033	BENCH: 600D, stainless steel lipped	①	1	[FIJO-021] with wet edge and integral sink
JOBE-231	BENCH SUPPORT: legs and frame, stainless steel	①	1	[FIJO-100]



Fittings, Furniture and Equipment (FF&E)

CODE	DESCRIPTION	GROUP	QTY	COMMENT
FIBM-019	WHITEBOARD: fixed, magnetic	②	1	[FIBM-019] optional
FIDI-106	DISPENSER: PPE, disposable gloves, wall mounted	②	1	[FIDI-014]
FIDI-231	DISPENSER: paper towel, wall mounted	②	1	[FIDI-025] to basin
FIDI-256	DISPENSER: soap, wall mounted	②	1	[FIDI-030] to basin
FIDI-261	DISPENSER: soap, antimicrobial, wall mounted	②	2	[FIDI-031] to basin
FIHR-029	HOOK: equipment, wall mounted	①	1	[FIHR-029] optional, for lifter straps/slings
FIHR-111	BRACKET: sharps bin, wall mounted	②	2	[FIHR-016]
FIHR-411	BRACKET: patient monitor, articulated, wall mounted	②	1	[FIHR-067]
FIRT-101	GRAB RAIL: straight, horizontal	①	2	[FIRT-005] to boot wash
FIRT-116	GRAB RAIL: U-shaped, floor mounted	①	1	[FIRT-046] to boot wash
FQBS-202	STOOL: clinical areas, mobile	③	1	[FQBS-051] optional, dependent on project requirements
FQCL-211	LINEN CARRIER: dirty, single skip	③	1	[MMGE-048]
FQGE-121	CLOCK: analogue display, battery operated	③	1	[FQGE-005]
FQWS-003	BIN: sharps, clinical, medium	③	1	[FQCW-043]
FQWS-003	BIN: sharps, clinical, medium	③	1	[FQCW-044]
FQWS-051	BIN: general waste, 20L	③	1	[FQCW-016] to basin
FQWS-057	BIN: general waste, 120L	③	1	[FQCW-015]
FQWS-237	BIN: clinical waste, 120L	③	1	[FQCW-004]
HYBA-101	BASIN: type A, handwashing	①	1	[HYBA-009]
HYBA-301	SINK: single bowl	①	1	[HYBA-055]
HYGE-031	TROUGH: boot wash, floor recessed	①	1	[HYGE-026]

CODE	DESCRIPTION	GROUP	QTY	COMMENT
HYGE-151	HOSE REEL: retractable	①	1	[HYTP-003]
HYSH-111	STATION: emergency, safety, shower and eye wash, wall mounted	①	1	[HYTP-023]
HYTP-062	TAPWARE: basin, tap set, fixture mounted, lever handles	①	2	[HYTP-038]
HYTP-176	TAPWARE: sink, tap set, fixture mounted, laboratory	①	1	[HYTP-034] to dissection table
HYTP-406	TAPWARE: spout, wall mounted	①	1	[HYTP-028] to boot wash
HYTP-416	TAPWARE: tap, wall mounted	①	1	[HYTP-030] to boot wash
HYTP-421	TAPWARE: hosecock, wall mounted	①	1	[HYTP-012] to hose reel
ITNE-076	TELEPHONE: handset, wall mounted	②	1	[ITNE-053] suitable for use in wet zone
ITNE-121	COMPUTER: single display screen, all-in-one	③	1	[ITNE-066] suitable for use in wet zone
ITNE-201	INPUT DEVICE: keyboard and mouse	③	1	[ITNE-067] suitable for use in wet zone
LIFX-131	LIGHT: procedure, single articulated arm, 1 light head, ceiling mounted	①	1	[LIFX-014] hardwired, a second light may be provided dependent on service requirements
MMBE-211	TRACK: patient lifter, fixed, ceiling mounted	①	1	[MMBE-071] optional, mobile lifter may be used, bariatric capacity to be confirmed based on service requirements
MMBE-251	CASSETTE MOTOR: patient lifter, for ceiling mounted system	①	1	[MMBE-071] optional, mobile lifter may be used, bariatric capacity to be confirmed based on service requirements
MMBE-266	SLING: patient lifter, supine	③	1	[MMBE-071] optional, mobile lifter may be used, bariatric capacity to be confirmed based on service requirements
MMGE-201	SCALES: digital, benchtop	③	1	[MMGE-185]
MMMT-001	TABLE: evisceration, autopsy, fixed, height adjustable	③	1	[MMLA-086] bariatric capacity to be determined based on service requirements, services on table to be considered: suction capabilities (built-in or provided), power to table, power outlets (if not on table), down draft (optional, built-in to table), water connection (built-in to table)
MMMT-006	TABLE: dissection, fixed, height adjustable	③	1	[MMLA-087] with integral sink
MMTR-071	TROLLEY: equipment, 490W x 490D x 900H nom	③	2	[MMGE-119]



Engineering Services

CODE	DESCRIPTION	GROUP	QTY	COMMENT
ELGP-131	GPO: single, emergency power, wall mounted	①	1	[ELGP-123] to computer. IP rated to suit AS/NZS wet zone requirements
ELGP-191	GPO: single, special, for equipment, wall mounted	①	8	[ELGP-106] 4x outlets on evisceration and dissection tables, dependent on service requirements and vendor specifications. Required to be emergency power/generator backed. IP rated to suit AS/NZS wet zone requirements
ELGP-191	GPO: single, special, for equipment, wall mounted	①	2	[ELGP-106] power to evisceration and dissection tables (if not part of table configuration). Required to be emergency power/generator backed. IP rated to suit AS/NZS wet zone requirements
ELGP-231	GPO: double, emergency power, wall mounted	①	2	[ELGP-223] IP rated to suit AS/NZS wet zone requirements
ELSW-001	SWITCH: light	①	2	[ELBO-015] 1x general room lights, 1x task light over dissection bench. IP rated to suit AS/NZS wet zone requirements
HYDR-021	DRAIN: floor waste, round	①	1	[HYDR-004] to boot wash
HYDR-025	DRAIN: floor waste, linear	①	2	[HYDR-003]

CODE	DESCRIPTION	GROUP	QTY	COMMENT
HYDR-211	DIRECT CONNECTION: wastewater, evisceration/autopsy table	①	2	[HYDR-001] 1 to autopsy table, 1 to evisceration table; confirm requirements at project level
HYTP-451	DIRECT CONNECTION: water, cold	①	2	[HYTP-008] connection only, to evisceration and dissection tables, dependent on service requirements and vendor specifications
HYTP-451	DIRECT CONNECTION: water, cold	①	2	1 to sink, 1 to basin
HYTP-451	DIRECT CONNECTION: water, cold	①	1	1 to hosecock, 1 to boot trough tapware
HYTP-452	DIRECT CONNECTION: water, hot	①	1	to sink
HYTP-452	DIRECT CONNECTION: water, hot	①	2	[HYTP-013] connection only, to evisceration and dissection tables, dependent on service requirements and vendor specifications
HYTP-453	DIRECT CONNECTION: water, warm	①	1	to basin
ITIN-016	OUTLET: data, single RJ45, wall mounted	①	3	[ITIN-016] IP rated to suit AS/NZS wet zone requirements
ITIN-026	OUTLET: data, double RJ45, wall mounted	①	1	[ITIN-026] IP rated to suit AS/NZS wet zone requirements
ITSE-166	INTERCOM: audio, slave, wall mounted	①	1	[ITSE-029] suitable for use in wet zone
LIFX-101	LIGHT: task, ceiling recessed	①	1	[LIFX-037] task lighting over dissection table, articulated task light is also acceptable dependent on service requirements
MEEX-122	DIRECT CONNECTION: exhaust, downdraft	①	1	[MEEX-017] optional, on evisceration and dissection tables, dependent on service requirements and vendor specifications
MGAS-063	OUTLET: suction, wall mounted	①	2	[MGAS-060] on evisceration and dissection tables, dependent on service requirements and vendor specifications

For guidance on how to use the Room Data Sheet, please visit the [AushFG website](#).