



Patient Bay - Medical Day Treatment

Room Code	PBTR-MD
Briefed Area	9.00 m ²
Ceiling Height	2.7 m
Occupancy	1 patient; 1 visitor; 1-2 staff intermittently
Hours of Operation	8am to 6pm, daily

Description | The Patient Bay - Medical Day Treatment may be used for assessment and treatment of non-admitted patients receiving scheduled treatment in either a bed or recliner (e.g. for infusions such as antibiotics and chemotherapy).

Electrical | PROTECTION: body protected
PROTECTION: cardiac protected

Lighting | LIGHTING: general
LIGHTING: colour corrected
LIGHTING: dimmable
LIGHTING: indirect

Nurse Call and Duress | NURSE CALL SYSTEM: buttons / handset

HVAC | AIRCONDITIONING: general
AIRCONDITIONING: HEPA filtered
AIRCONDITIONING: positive pressure
AIRCONDITIONING: negative pressure
VENTILATION: exhaust
VENTILATION: natural

Additional Considerations

- A staff handwashing basin will be provided for each 4 bays
- The Patient Bay will have, wherever possible, clear line of sight to a staff station and close access to a patient toilet.
- Access to natural light and views is required.
- Medical Gases may be provided individually to each bay or shared between two bays depending on service requirements.



Doors and Windows

CODE	DESCRIPTION	COMMENT
WIFX-162.01	WINDOW: fixed, external, double glazed, sill at 1650H	[DWWH-005] access to natural light is required
AHWCV-001	WINDOW COVERING: blind, roller, single	[DWWD-004] dependent on local infection control policies and service requirements



Finishes

CODE	DESCRIPTION	COMMENT
FLVY-101	FLOOR FINISH: vinyl, seamless, standard slip resistance	[FLVY-007]
FLSK-021	SKIRTING: vinyl, integral with floor vinyl, covered	[FLSK-010]
WLFI-001	WALL FINISH: paint	[WLWA-004]
CLTI-022.02	CEILING: drop-in tiles, acoustic, prefinished, 600 x 1200	[CLTI-002]
CLCN-041	CORNICE: wall trim fixing, prefinished	[CLCN-009]
WLPR-002	WALL PROTECTION: bedhead wall panel, low height	[WLPR-002]



Fittings, Furniture and Equipment (FF&E)

CODE	DESCRIPTION	GROUP	QTY	COMMENT
FIDI-241	DISPENSER: alcohol-based hand rub, wall mounted	2	1	[FIDI-001]
FIHR-101	BRACKET: suction bottle, wall mounted	2	1	[FIHR-019] may be shared between two bays depending on service requirements
FIHR-111	BRACKET: sharps bin, wall mounted	2	1	[FIHR-016]
FIHR-214	BRACKET: computer display screen, single, articulated, ceiling mounted	2	1	[FIHR-004]
FIRT-051	TRACK: curtain, privacy	1	1	[FIRT-032]
FQBS-101	CHAIR: visitor, clinical areas	3	1	[FQBS-030]
FQBS-202	STOOL: clinical areas, mobile	3	1	[FQBS-051] quantity to be rationalised across the unit to suit number of staff
FQGE-101	CURTAIN: privacy screen	3	2	[FQGE-010]
FQSN-101	BEDSIDE TABLE: floor standing, patient locker	3	1	[FQSN-051]
FQTA-091	TABLE: overbed, mobile, height adjustable	3	1	[FQTA-013]
FQWS-003	BIN: sharps, clinical, medium	3	1	[FQCW-043]
ITCL-301	DISPLAY SCREEN: patient entertainment system	2	1	[ITBU-012]
LIFX-113	LIGHT: examination, single articulated arm, 1 light head, ceiling mounted	1	1	[LIFX-027] optional, dependent on service requirements
MGFP-061	ADAPTER: suction	3	1	[MMGE-002] may be shared between two bays depending on service requirements
MMBE-091	CHAIR: patient, recliner, clinical, electric	3	1	[MMGE-023] or bed
MMGE-191	CANNISTER: suction bottle	3	1	[MMGE-100] may be shared between two bays depending on service requirements
MMTR-071	TROLLEY: equipment, 490W x 490D x 900H nom	3	1	[MMGE-119]



Engineering Services

CODE	DESCRIPTION	GROUP	QTY	COMMENT
ELGP-101	GPO: single, wall mounted	1	1	[ELGP-109] low level, for bed/recliner
ELGP-106	GPO: single, ceiling mounted	1	1	[ELGP-101] to patient entertainment screen
ELSW-001	SWITCH: light	1	1	[ELBO-015] to exam light
ITCL-191	BUTTON: nurse call, emergency, with cancel, wall mounted	1	1	[ITCL-004]
ITIN-071	OUTLET: MATV, ceiling mounted	1	1	[ITIN-070] to patient entertainment screen, may be replaced with data outlet in accordance with ICT policy
LIFX-161	LIGHT: reading, wall mounted	1	1	[LIFX-011] hardwired



Medical Services Panels and Pendants

CODE	DESCRIPTION	GROUP	QTY	COMMENT
Panel 1				
MMSP-051	MEDICAL SERVICES PANEL: wall mounted	1	1	[MMSP-026]
ELGP-216	GPO: double, on services panel	1	3	[ELGP-206]
ELPR-083	RCD: residual current device, emergency power, on services panel	1	1	[ELPR-006]

CODE	DESCRIPTION	GROUP	QTY	COMMENT
ITCL-133	BUTTON: nurse call, patient to staff call, with cancel, with handset connection point, on services panel	①	1	[ITCL-006] + indicator button & light
ITCL-183	BUTTON: nurse call, staff assist, with cancel, on services panel	①	1	[ITCL-008] + indicator button & light
ITCL-252	HANDSET: nurse call, patient to staff call, with light control	①	1	[ITCL-028] call button on the handset, reading light control
ITIN-025	OUTLET: data, double RJ45, on services panel	①	1	[ITIN-025]
Panel 2				
MMSP-051	MEDICAL SERVICES PANEL: wall mounted	①	1	[MMSP-026] 1 panel for medical gases may be shared between two bays depending on service requirements
MGAS-022	OUTLET: medical air, on services panel	①	1	[MGAS-022] may be shared between two bays depending on service requirements
MGAS-042	OUTLET: oxygen (O2), on services panel	①	1	[MGAS-042] may be shared between two bays depending on service requirements
MGAS-062	OUTLET: suction, on services panel	①	1	[MGAS-062] may be shared between two bays depending on service requirements
MGFP-021	FLOWMETER: medical air	③	1	[MGFP-002] may be shared between two bays depending on service requirements
MGFP-041	FLOWMETER: oxygen	③	1	[MGFP-004] may be shared between two bays depending on service requirements
MGFP-061	ADAPTER: suction	③	1	[MMGE-002] may be shared between two bays depending on service requirements

For guidance on how to use the Room Data Sheet, please visit the [AusHFG website](#).