



Plaster Room

Room Code	PLST
Briefed Area	18.00 m ²
Ceiling Height	2.7 m
Occupancy	1 patient; 1 visitor; up to 4 staff
Hours of Operation	Up to 24 hours (dependent on unit operating hours)

Description | The Plaster Room allows for the application of plaster or other splints and for the closed reduction of displaced fractures or dislocations under sedative or regional anaesthesia.

Electrical | PROTECTION: body protected
PROTECTION: cardiac protected

Lighting | LIGHTING: general
LIGHTING: colour corrected
LIGHTING: dimmable
LIGHTING: indirect

Nurse Call and Duress | NURSE CALL SYSTEM: buttons / handset

HVAC | AIRCONDITIONING: general
AIRCONDITIONING: HEPA filtered
AIRCONDITIONING: positive pressure
AIRCONDITIONING: negative pressure
VENTILATION: exhaust
VENTILATION: natural

Additional Considerations

- Typically provided as a dedicated room for high volume services, otherwise, a generic Treatment Room may be used with the inclusion of a plaster trap and other specific plastering equipment.
- The cast saw shall incorporate an integrated dust extraction system or the room shall have an independent dust extraction system suitably located to capture the dust at source. A room exhaust shall also be included to prevent dust from entering adjacent spaces.
- The plaster trap shall be easily accessible for cleaning and maintenance.



Doors and Windows

CODE	DESCRIPTION	COMMENT
DOHI-114.01	DOOR: hinged, 1 1/2 leaves, 1400 clear opening, solid, standard vision panel	[DWGL-033]
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to corridor side, extent dependent on movement of beds and mobile equipment in adjacent area; to be coordinated with corridor wall protection
AFDPR-006.01	DOOR PROTECTION: plate, to 900H	[DWPR-005] to room side
AFDPR-056.03	DOOR FRAME PROTECTION: full wrap, to 1200H	[DWPR-010] optional, provision and extent dependent on frame material/finish and movement of beds and mobile equipment in adjacent area



Finishes

CODE	DESCRIPTION	COMMENT
FLVY-101	FLOOR FINISH: vinyl, seamless, standard slip resistance	[FLVY-007] floor finish must be minimum slip rating R10 / Pendulum P3 or agreed equivalent to extend under basin or inset section of non-slip vinyl required

CODE	DESCRIPTION	COMMENT
FLSK-021	SKIRTING: vinyl, integral with floor vinyl, covered	[FLSK-010]
WLFI-001	WALL FINISH: paint	[WLWA-004]
WLFI-011.04	WALL FINISH: vinyl, to 1500 AFFL	[WLWA-020]
WLFI-101	SPLASHBACK: vinyl	[WLWA-007] to basin
CLFS-011	CEILING: flush set, suspended	[CLFS-006]
CLFI-002	CEILING FINISH: paint, clinical areas	[CLFS-006]
CLCN-031	CORNICE: square set	[CLCN-008]



Joinery

CODE	DESCRIPTION	GROUP	QTY	COMMENT
JOBE-033	BENCH: 600D, stainless steel lipped	①	1	[FIJO-060] with sink and wet edge
JOCU-121	CUPBOARD: wall mounted, double door	①	3	[FIJO-145]
JOCU-221	CUPBOARD: under bench, double door	①	3	[FIJO-136]
JOGE-001	BULKHEAD: joinery	①	1	[FIJO-212] plasterboard bulkhead also acceptable



Fittings, Furniture and Equipment (FF&E)

CODE	DESCRIPTION	GROUP	QTY	COMMENT
FIDI-106	DISPENSER: PPE, disposable gloves, wall mounted	②	1	[FIDI-014]
FIDI-231	DISPENSER: paper towel, wall mounted	②	1	[FIDI-025] to basin
FIDI-241	DISPENSER: alcohol-based hand rub, wall mounted	②	1	[FIDI-001]
FIDI-256	DISPENSER: soap, wall mounted	②	1	[FIDI-030] to basin
FIDI-261	DISPENSER: soap, antimicrobial, wall mounted	②	1	[FIDI-031] to basin
FIHR-101	BRACKET: suction bottle, wall mounted	②	2	[FIHR-019] quantity to suit number of suction outlets
FIRT-051	TRACK: curtain, privacy	①	1	[FIRT-039]
FQGE-101	CURTAIN: privacy screen	③	1	[FQGE-010]
FQWS-051	BIN: general waste, 20L	③	1	[FQCW-016] to basin
FQWS-233	BIN: clinical waste, 50L	③	1	[FQCW-008] optional; type, size and quantity of bins to be determined based on service requirements and waste management policies
HYBA-121	BASIN: type B, handwashing	①	1	[HYBA-010]
HYBA-301	SINK: single bowl	①	1	[HYBA-053] hot water
HYTP-062	TAPWARE: basin, tap set, fixture mounted, lever handles	①	1	[HYTP-038] to sink
HYTP-067	TAPWARE: basin, tap set, wall mounted, lever handles	①	1	[HYTP-037] to basin
ITBU-441	CLOCK: analogue display, electric, wall mounted	①	1	[ITCL-014]
ITNE-076	TELEPHONE: handset, wall mounted	②	1	[ITNE-053]
ITNE-111	COMPUTER: workstation on wheels (WOW), medical grade, all-in-one computer, keyboard and mouse, on cart	③	1	[ITNE-008+MMGE-164] computer may alternatively be on benchtop or wall mounted; allocation of WOWs across department to be rationalised
LIFX-113	LIGHT: examination, single articulated arm, 1 light head, ceiling mounted	①	1	[LIFX-027]

CODE	DESCRIPTION	GROUP	QTY	COMMENT
MMBE-156	TROLLEY: patient, examination	③	1	[MMBE-043] or plinth/examination couch
MMGE-191	CANNISTER: suction bottle	③	2	[MMGE-100] quantity to suit number of suction outlets
MMGE-233	DIAGNOSTIC SET: otoscope and ophthalmoscope, wall mounted	②	1	[MMGE-031] optional, may be portable/mobile on stand
MMGE-326	DEVICE: saw, plaster, with integral dust extraction	③	1	[MMAH-023 + MMAH-045] integral/localised dust extraction required or independent dust extraction system to be provided
MMTR-141	TROLLEY: plaster	③	1	[MMAH-028]



Engineering Services

CODE	DESCRIPTION	GROUP	QTY	COMMENT
ELGP-101	GPO: single, wall mounted	①	1	[ELGP-109] low level, for patient trolley
ELGP-105	GPO: single, wall mounted, cleaner	①	1	[ELGP-110] as required for distribution across the unit in accordance with AS/NZS 3003
ELGP-106	GPO: single, ceiling mounted	①	1	[ELGP-101] to exam light
ELGP-201	GPO: double, wall mounted	①	4	[ELGP-208]
ELPR-071	RCD: residual current device, wall mounted	①	4	[ELPR-006]
ELSW-001	SWITCH: light	①	3	[ELBO-015] 1 x general room lighting, 1 x task light, 1 x examination light
HYDR-151	TRAP: plaster	①	1	[HYDR-008] to be easily accessible for cleaning and maintenance
HYDR-201	DIRECT CONNECTION: wastewater, general	①	2	1 to sink, 1 to basin
HYTP-451	DIRECT CONNECTION: water, cold	①	2	1 to sink tapware, 1 to basin tapware
HYTP-452	DIRECT CONNECTION: water, hot	①	1	to sink tapware
HYTP-453	DIRECT CONNECTION: water, warm	①	1	to basin tapware
ITCL-111	LIGHT: nurse call indicator, ceiling mounted	①	1	outside entry door
ITCL-191	BUTTON: nurse call, emergency, with cancel, wall mounted	①	1	[ITCL-004] + indicator button & light
ITCL-251	HANDSET: nurse call, patient to staff call	①	1	[ITCL-028] call button on the handset
ITIN-016	OUTLET: data, single RJ45, wall mounted	①	1	[ITIN-016] to telephone
ITIN-026	OUTLET: data, double RJ45, wall mounted	①	1	[ITIN-026]
ITSE-061	BUTTON: security, duress, fixed, wall mounted	①	1	[ITSE-006]
LIFX-108	LIGHT: task, integral to joinery	①	3	[LIFX-008]



Medical Services Panels and Pendants

CODE	DESCRIPTION	GROUP	QTY	COMMENT
Panel 1				
MMSP-051	MEDICAL SERVICES PANEL: wall mounted	①	1	[MMSP-026]
ELGP-116	GPO: single, on services panel	①	3	[ELGP-105]
ELGP-146	GPO: single, emergency power, on services panel	①	3	[ELGP-122]
ITCL-133	BUTTON: nurse call, patient to staff call, with cancel, with handset connection point, on services panel	①	1	[ITCL-006] + indicator button & light
ITCL-183	BUTTON: nurse call, staff assist, with cancel, on services panel	①	1	[ITCL-008] + indicator button & light

CODE	DESCRIPTION	GROUP	QTY	COMMENT
ITIN-025	OUTLET: data, double RJ45, on services panel	①	1	[ITIN-025]
MGAS-022	OUTLET: medical air, on services panel	①	1	[MGAS-022]
MGAS-032	OUTLET: nitrous oxide (N2O), on services panel	①	1	[MGAS-032] optional
MGAS-042	OUTLET: oxygen (O2), on services panel	①	1	[MGAS-042] 1 is optional, quantity dependent on service requirements
MGAS-042	OUTLET: oxygen (O2), on services panel	①	1	[MGAS-042]
MGAS-052	OUTLET: scavenge, on services panel	①	1	[MGAS-053]
MGAS-062	OUTLET: suction, on services panel	①	1	[MGAS-062] 1 is optional, quantity dependent on service requirements
MGAS-062	OUTLET: suction, on services panel	①	1	[MGAS-062]
MGFP-021	FLOWMETER: medical air	③	1	[MGFP-002]
MGFP-041	FLOWMETER: oxygen	③	2	[MGFP-004] quantity to suit number of oxygen outlets
MGFP-051	ADAPTER: scavenge	③	1	[MMGE-001] optional
MGFP-061	ADAPTER: suction	③	1	[MMGE-002] quantity to suit number of suction outlets
MGFP-061	ADAPTER: suction	③	1	[MMGE-002] optional; quantity to suit number of suction outlets
MGFP-082	BLENDER: oxygen (O2) and nitrous oxide (N2O)	③	1	[MMMA-011] optional, may be mobile, on trolley or wall mounted

For guidance on how to use the Room Data Sheet, please visit the [AushFG website](#).