



Room Name	Clean-Up Room - Shared, 15m2	Room Number CLUP-15			
Revision	4				
Issue Date	20.03.17				
BRIEFING					
Briefed Area	15.00 m ²				
Hours of Operation	24 Hours				
Occupancy	1 - 2				
Acoustics	· · -				
Remarks					
Description / Special Requirements	The Clean-Up Room, Shared is used for rinsing and cleaning of medical equipment or endoscopes used in the Operating Rooms or Procedure Rooms. Low temperature sterilisation of endoscopes may also be undertaken in this room. The shared Clean-up Room should be located between two Minor Operating Rooms or Procedures rooms where endoscopy procedures are undertaken.				
SERVICES					
General Requirements	BODY PROTECTED = No	CARDIAC PROTECTED = No			
Air	AIRCONDITIONING = Yes	AIRCONDITIONING: hepa filtered = No			
	AIRCONDITIONING: positive pressure = No	AIRCONDITIONING: negative pressure = No			
	EXHAUST: room exhaust = Yes	VENTILATION = No			
Lighting	LIGHTING: colour corrected = No	LIGHTING: dimmable = No			
Ligiting	LIGHTING: downlights = No	LIGHTING: downlights, dimmable = No			
	LIGHTING: fluorescent, downlights = No	LIGHTING: fluorescent, flush with ceiling, tamper proof = No			
	LIGHTING: fluorescent, general = Yes	LIGHTING: indirect = No			
Nurse Call	NURSE CALL SYSTEM = No				

Room Fabrics

AusHFG code	Description	Comment
CLTI-014	CEILING: plasterboard, water resistant, drop-in tiles, prefinished, 600 x 1200	
CLCN-002	CORNICE: aluminium, 24mm shadow angle, powdercoat	
DWGL-020	DOOR: 1100mm c/o, single leaf, observation panel, solid core/glass, paint	1100 mm clear opening to ORs, double swing, sliding cover to observation panel
DWGL-005	DOOR: 910mm c/o, single leaf, observation panel, solid core/glass, paint	
DWPR-006	DOOR PROTECTION: protection plate, to 900 AFFL, PVC, prefinished	Equal to Acrovyn
DWPR-010	DOOR PROTECTION: to door frame, PVC, prefinished	Equal to Acrovyn
FLVY-004	FLOOR FINISH: vinyl, seamless, coved, non-slip, safety	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLWA-020	WALL FINISH: vinyl, to 1500 AFFL	

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	1	for tube cleaning
HYBA-010	BASIN: handwash, Type B	1	1	
FIJO-032	BENCH: 750D, open under, stainless steel lipped	1	1	with wet edge and inset sinks

AHHA

australisian health infrastructur alliance

distributed in automatical in structure and infrastructure and infrastructure

Page: 1 Print Date: 20/03/2017





Room Name

Clean-Up Room - Shared, 15m2

Room Number CLUP-15

AusHFG code	Description	Group	Qty	Comment
FIJO-059	BENCH TOP: 600D, joinery under, laminate finish	1	1	
FQCW-006	BIN: waste, clinical, 240L	3	1	
FQCW-016	BIN: waste, general, 20L	3	1	
FIHR-016	BRACKET: for sharps bin	2	1	
FIJO-104	CUPBOARD: drawer unit, 3 equal	1	3	under dry bench
FIJO-143	CUPBOARD: wall mount, double door, adj shelves	1	2	over dry bench
FIJO-156	CUPBOARD: wall mount, single door, adj shelves	1	1	over dry bench
FIDI-014	DISPENSER: disposable glove	2	2	
FIDI-025	DISPENSER: paper/hand towel	2	1	to basin
FIDI-030	DISPENSER: soap	2	1	to basin
FIDI-031	DISPENSER: soap antiseptic	2	1	to basin
MGFP-002	FLOWMETER: medical air	3	1	for tube drying
HYTP-006	HOSE: spray, pressure control	1	1	to flushing rim sink, with flush button
FIJO-196	SHELF: stainless steel, overbench, wall mount	1	1	over sink bench
FIJO-197	SHELF: stainless steel, underbench	1	1	under sink bench
HYBA-046	SINK: flushing rim	1	1	in bench
HYBA-063	SINK: single, 750 x 500 x 300D	1	1	hot water
MMIC-054	STERILISER: low temperature, rapid cycle, single door	3	1	optional
HYTP-037	TAPWARE: mixer, elbow levers	1	2	1 to sink, 1 to basin

Engineering Services

Engineering der vices				
AusHFG code	Description	Group	Qty	Comment
HYDR-004	DRAIN: floor waste	1	1	
HYDR-006	DRAIN: tundish	1	1	for low temp steriliser, optional
ELGP-208	GPO: double, wall mounted	1	1	
ELGP-223	GPO: emergency power, double, wall mounted	1	1	
ELGP-106	GPO: single, special, for equipment	1	1	optional, for low temp steriliser
LIFX-008	LIGHT: task light, built in	1	3	under O/H cupboards
HYTP-008	OUTLET: cold water	1	1	for low temp steriliser, optional
HYTP-013	OUTLET: hot water	1	1	for low temp steriliser, optional
MGAS-023	OUTLET: medical air (MA), wall mounted	1	1	over dry bench for tube drying
MGAS-063	OUTLET: suction, wall mounted	1	1	near sink for tube cleaning
ELBO-015	SWITCH: light, single	1	3	

Page: 2 Print Date: 20/03/2017