

Room Name	Patient Bay, Emergency - Acute Tr	eatm	ent Room Number PBTR-A	
Revision	5			
Issue Date	29.05.2019			
BRIEFING				
Briefed Area	40.00 2			
	12.00 m <sup>2</sup> 24 Hours			
Hours of Operation Occupancy	24 Hours 1 patient; 1-2 staff			
Acoustics	i patient, 1-2 stan			
Remarks				
Description / Special Requirements	The Patient Bay, Emergency - Acute Trea	າ the E	t may be used for assessment and treatment of patients w Emergency Unit. The bay may have partition side nib walls	ith . The
Opecial Requirements		aff ha	ndwashing basin and, wherever possible, clear line of sigh	
SERVICES	Patient Bay will have close access to a st	aff ha	ndwashing basin and, wherever possible, clear line of sigh	
	Patient Bay will have close access to a st	aff ha	ndwashing basin and, wherever possible, clear line of sigh  CARDIAC PROTECTED	t to a
SERVICES	Patient Bay will have close access to a st staff station.			t to a
SERVICES General Requirements	Patient Bay will have close access to a st staff station.  BODY PROTECTED	<b>√</b>	CARDIAC PROTECTED	t to a
SERVICES General Requirements	Patient Bay will have close access to a st staff station.  BODY PROTECTED  AIRCONDITIONING	<b>V</b>	CARDIAC PROTECTED  AIRCONDITIONING: hepa filtered	t to a
SERVICES General Requirements	Patient Bay will have close access to a st staff station.  BODY PROTECTED  AIRCONDITIONING  AIRCONDITIONING: positive pressure	<b>V</b>	CARDIAC PROTECTED  AIRCONDITIONING: hepa filtered  AIRCONDITIONING: negative pressure	
SERVICES General Requirements Air	Patient Bay will have close access to a st staff station.  BODY PROTECTED  AIRCONDITIONING  AIRCONDITIONING: positive pressure EXHAUST: room exhaust	<b>V</b>	CARDIAC PROTECTED  AIRCONDITIONING: hepa filtered  AIRCONDITIONING: negative pressure  VENTILATION	t to a
SERVICES General Requirements Air	Patient Bay will have close access to a st staff station.  BODY PROTECTED  AIRCONDITIONING AIRCONDITIONING: positive pressure EXHAUST: room exhaust  LIGHTING: colour corrected	✓ ✓ □	CARDIAC PROTECTED  AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure VENTILATION  LIGHTING: dimmable LIGHTING: downlights, dimmable LIGHTING: fluorescent/LED, flush with ceiling, tamper	t to a
SERVICES General Requirements Air	Patient Bay will have close access to a st staff station.  BODY PROTECTED  AIRCONDITIONING AIRCONDITIONING: positive pressure EXHAUST: room exhaust  LIGHTING: colour corrected LIGHTING: downlights	<b>V</b>	CARDIAC PROTECTED  AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure VENTILATION  LIGHTING: dimmable LIGHTING: downlights, dimmable	t to a

## **Room Fabrics**

AusHFG code	Description	Comment
CLFS-006	CEILING: plasterboard, flush set, suspended, paint, washable	
CLCN-008	CORNICE: square set flush between wall and ceiling, paint	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLWA-025	WALL FINISH: vinyl, to 1200 AFFL	
WLPR-002	WALL PROTECTION: bed protection wall panel, composite, laminate finish	high pressure laminate concealed fixing
WLPR-013	WALL PROTECTION: corner guards, to 1500 AFFL, prefinished PVC	as required

## Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	2	
MMBE-011	BED: patient, emergency, electric	3	1	
FQSN-006	BIN: storage bins, for wall mounted panel	3	3	wire basket bins
FQCW-005	BIN: waste, clinical, 20L	3	1	
FQCW-016	BIN: waste, general, 20L	3	1	

AHHA auditation booth infrastructure allicore

Page: 1 Print Date: 27/05/2019





**Room Name** 

**Patient Bay, Emergency - Acute Treatment** 

**Room Number PBTR-A** 

AusHFG code	Description	Group	Qty	Comment
FQCW-044	BIN: waste, sharps, wall mounted	3	1	
FIHR-081	BRACKET: for all-in-one computer, with keyboard shelf, wall mounted, articulated	2	1	optional
FIHR-067	BRACKET: for monitor, wall mounted, articulated, adj height	2	1	
FIHR-016	BRACKET: for sharps bin	2	1	
FIHR-019	BRACKET: for suction bottle	2	2	
FQBS-030	CHAIR: visitor, patient	3	1	
ITNE-066	COMPUTER: all-in-one, wall mounted	3	1	optional, or workstation on wheels
ITNE-008	COMPUTER: for workstation on wheels	3	1	
FQGE-009	CURTAIN: bed screen	3	2	
MMGE-031	DIAGNOSTIC SET: wall mounted	2	1	
FIDI-001	DISPENSER: alcohol based hand rub	2	1	
FIDI-014	DISPENSER: disposable glove	2	1	
MGFP-002	FLOWMETER: medical air	3	1	
MGFP-004	FLOWMETER: oxygen	3	2	
ITNE-067	KEYBOARD: keyboard and mouse only	3	1	optional
LIFX-013	LIGHT: procedure, ceiling mounted, small	1	1	
MMHA-005	MONITOR: patient, cardiac	3	1	
FIHR-033	PANEL: for storage bins, wall mount	1	1	
MMGE-100	SUCTION BOTTLE	3	2	
FQTA-013	TABLE: overbed	3	1	
FIRT-025	TRACK: curtain, bedscreen, U shaped	1	1	
MMGE-164	WORKSTATION ON WHEELS: clinical, mobile	3	1	

## **Engineering Services**

AusHFG code	Description	Group	Qty	Comment
ITCL-004	BUTTON: nurse call, emergency	1	1	on services panel, + indicator button & light
ITCL-006	BUTTON: nurse call, patient to staff	1	1	on services panel, + indicator button & light
ITCL-008	BUTTON: nurse call, staff assist	1	1	on services panel, + indicator button & light
ELGP-223	GPO: emergency power, double, wall mounted	1	1	for WOW/ computer
ELGP-125	GPO: emergency power, single, ceiling mounted	1	1	
ELGP-122	GPO: emergency power, single, on services panel	1	8	on services panel
ELGP-123	GPO: emergency power, single, wall mounted	1	3	optional, for cardiac monitor and wall mounted PC, if required
ELGP-110	GPO: single, wall mounted, cleaner	1	1	shared between a group of beds
ITCL-028	HANDSET: nurse call, with patient to staff call button	1	1	call button on the handset
ITIN-025	OUTLET: data, double RJ45, on services panel	1	1	

AHHA

Page: 2 Print Date: 27/05/2019



**RDS Templates** 



**Room Name** 

**Patient Bay, Emergency - Acute Treatment** 

**Room Number PBTR-A** 

AusHFG code	Description	Group	Qty	Comment
ITIN-026	OUTLET: data, double RJ45, wall mounted	1	2	optional, for cardiac monitor and wall mounted PC, if required
MGAS-022	OUTLET: medical air (MA), on services panel	1	1	on services panel
MGAS-042	OUTLET: oxygen (O2), on services panel	1	2	on services panel
MGAS-062	OUTLET: suction, on services panel	1	2	on services panel
ELPR-006	RCD: residual current device	1	2	
MMSP-026	SERVICES PANEL: wall mounted	1	3	
ELBO-015	SWITCH: light, single	1	1	

