



Room Name	SPECT-CT Control Room		Room Number SPECT-CTCR	
Revision	1			
Issue Date	29.05.2019			
BRIEFING				
Briefed Area	14.00 m ²			
Hours of Operation	Available 24 hours			
Occupancy	Up to 3 staff			
Acoustics	Acoustic privacy required			
Remarks				
Description / Special Requirements	Room adjacent to the SPECT-CT Imaging equipment, and review of procedure image	g Roo ges. A	m with direct visibility to the patient for remote control of dditional Design Considerations: X-Ray shielding as advise	ed by
	Consultant. Cable ducts from the imaging	room	and equipment room will be required.	,
	Consultant. Cable ducts from the imaging	room	and equipment room will be required.	
SERVICES General Requirements	Consultant. Cable ducts from the imaging	room	and equipment room will be required. CARDIAC PROTECTED	
SERVICES	Consultant. Cable ducts from the imaging	room	CARDIAC PROTECTED	
SERVICES General Requirements	Consultant. Cable ducts from the imaging BODY PROTECTED AIRCONDITIONING	room	CARDIAC PROTECTED AIRCONDITIONING: hepa filtered	
SERVICES General Requirements	Consultant. Cable ducts from the imaging	y room	CARDIAC PROTECTED	
SERVICES General Requirements	Consultant. Cable ducts from the imaging BODY PROTECTED AIRCONDITIONING AIRCONDITIONING: positive pressure	y room	CARDIAC PROTECTED AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure	
SERVICES General Requirements Air	BODY PROTECTED AIRCONDITIONING AIRCONDITIONING: positive pressure EXHAUST: room exhaust	y room	CARDIAC PROTECTED AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure VENTILATION	
SERVICES General Requirements Air	BODY PROTECTED AIRCONDITIONING AIRCONDITIONING: positive pressure EXHAUST: room exhaust LIGHTING: colour corrected	v .	CARDIAC PROTECTED AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure VENTILATION LIGHTING: dimmable	
SERVICES General Requirements Air	BODY PROTECTED AIRCONDITIONING AIRCONDITIONING: positive pressure EXHAUST: room exhaust LIGHTING: colour corrected LIGHTING: downlights	y com	CARDIAC PROTECTED AIRCONDITIONING: hepa filtered AIRCONDITIONING: negative pressure VENTILATION LIGHTING: dimmable LIGHTING: downlights, dimmable LIGHTING: fluorescent/LED, flush with ceiling, tamper	

Room Fabrics

AusHFG code	Description	Comment
CLTI-002	CEILING: acoustic, drop-in tiles, 600 x 1200	
CLCN-002	CORNICE: aluminium, 24mm shadow angle, powdercoat	
DWGL-005	DOOR: 910mm c/o, single leaf, observation panel, solid core/glass, paint	
DWSC-005	DOOR: 910mm c/o, single leaf, solid core, lead shielding, paint	
DWPR-005	DOOR PROTECTION: protection plate to 900 AFFL	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
DWOB-009	OBSERVATION: internal window, lead glass, clear	safety glass
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
FIJO-031	BENCH: 750D, open under, laminate finish	1	2	high bench at 900-950mm AFFL optional
FQCW-016	BIN: waste, general, 20L	3	1	
DWWD-004	BLIND: roller	3	1	to observation window for patient privacy

AHHA
autrition hath internation allows
conserve transports to the property tof the property to the property to the property to the property to

Page: 1 Print Date: 27/05/2019





Room Name

SPECT-CT Control Room

Room Number SPECT-CTCR

AusHFG code	Description	Group	Qty	Comment
FIHR-075	BRACKET: for computer CPU, underdesk mounted	2	2	
FIHR-076	BRACKET: for display screen, triple, desk mounted	3	1	
FIHR-066	BRACKET: for display screen, wall mounted, articulated, adj height	2	1	
FQSN-031	CABINET: short, single door, tambour, adj shelves, lockable	3	1	
FIJO-099	CAC: cable access cap	3	2	
FQBS-018	CHAIR: office, ergonomic	3	2	
ITNE-068	COMPUTER: hi performance, for clinical / PACS use, CPU, underbench mount	3	1	
MMMI-269	COMPUTER: tower, acquisition, SPECT/CT	1	1	
MMMI-274	CONSOLE: control, SPECT/CT	1	1	
FIDI-001	DISPENSER: alcohol based hand rub	2	1	
MMMI-369	DISPLAY SCREEN: contrast injector	3	1	
ITNE-063	DISPLAY SCREEN: for computer	3	1	
ITNE-069	DISPLAY SCREEN: PACS, portrait	3	2	
MMGE-032	DISPLAY SCREEN: slave, for patient monitor	3	1	optional
FIHR-026	HOOK: coat	1	2	optional
ITNE-067	KEYBOARD: keyboard and mouse only	3	1	
ITAV-014	MICROPHONE: surface mounted, gooseneck	2	1	communication with patient, on bench required if integrated intercom system not provided by equipment manufacturer
ITNE-023	PRINTER: networked	3	1	on benchtop
FIHR-044	RACK: lead apron	1	1	
ITAV-057	SPEAKER: ceiling recessed	2	1	communication with patient, on bench required if integrated intercom system not provided by equipment manufacturer
ITNE-049	TELEPHONE: handset, desktop	3	1	
FIBM-019	WHITEBOARD: fixed, magnetic	2	1	

Engineering Services

AusHFG code	Description	Group	Qty	Comment
ELBO-001	BUTTON: emergency stop	1	1	mains power, consult manufacturer for safety and compatibility with equipment.
ITCL-004	BUTTON: nurse call, emergency	1	1	
ELGP-208	GPO: double, wall mounted	1	2	one for printer
ELGP-223	GPO: emergency power, double, wall mounted	1	4	
ELGP-123	GPO: emergency power, single, wall mounted	1	1	
ELGP-110	GPO: single, wall mounted, cleaner	1	1	
ELGP-135	GPO: UPS, single, wall mounted	1	5	
ITSE-028	INTERCOM: audiovisual, wall mounted	1	1	to screening room bed entry door
LIFX-002	LIGHT: illuminated sign, "X-Ray in Use"	1	1	auto on when modality is in use



Page: 2 Print Date: 27/05/2019 Generated from dRofus © 2017 dRofus Pty Ltd



RDS Templates



Room Name SPECT-CT Control Room

Room Number SPECT-CTCR

AusHFG code	Description	Group	Qty	Comment
ITIN-026	OUTLET: data, double RJ45, wall mounted	1	4	
ITIN-016	OUTLET: data, single RJ45, wall mounted	1	3	1 outlet for display screen, patient monitor slave; one outlet for telephone
ELBO-012	SWITCH: light, dimmer	1	2	

