

Room Name	Uptake Room	Room Number UTRM
Revision	1	
Issue Date	24.06.2022	
BRIEFING		
Briefed Area	10.00 m ²	
Hours of Operation	8am to 6pm, daily	
Occupancy	1 patient, 1 staff	
Acoustics		
Remarks		
Description / Special Requirements	<p>The Uptake Room is used to inject patient with a radiopharmaceutical and allow them to rest until 'uptake' has occurred before transfer to the scanning room. Patients need to be relaxed in a warm, quiet and dimmed environment to promote uptake. Uptake duration is typically 45 to 60 minutes.</p> <p>Additional Considerations:</p> <ul style="list-style-type: none"> Visual and audio observation of the patient during this time is typically provided from the Control Room of the PET-CT Imaging Room and intercoms, cameras and display screens are to be provided to support patient observation. Injections of radiopharmaceuticals can be undertaken within the room or from outside the room via an automated dose delivery system, which will be accommodated outside each uptake room with a suitable lead-lined injection port in the wall that will allow the coil to be connected to both the machine and the patient (details should be confirmed with vendor) Radiation shielding is required to the room and is to be in accordance with a radiation shielding consultants advice. This will be based on the type of radiopharmaceutical, the number of patients seen on a daily basis and proximity to other services. 	
SERVICES		
Electrical	BODY PROTECTED <input checked="" type="checkbox"/>	CARDIAC PROTECTED <input type="checkbox"/>
HVAC	AIRCONDITIONING <input checked="" type="checkbox"/>	AIRCONDITIONING: hepa filtered <input type="checkbox"/>
	AIRCONDITIONING: positive pressure <input type="checkbox"/>	AIRCONDITIONING: negative pressure <input type="checkbox"/>
	EXHAUST: room exhaust <input type="checkbox"/>	VENTILATION <input type="checkbox"/>
Lighting	LIGHTING: colour corrected <input type="checkbox"/>	LIGHTING: dimmable <input checked="" type="checkbox"/>
	LIGHTING: downlights <input type="checkbox"/>	LIGHTING: downlights, dimmable <input type="checkbox"/>
	LIGHTING: fluorescent/LED, downlights <input type="checkbox"/>	LIGHTING: fluorescent/LED, flush with ceiling, tamper proof <input type="checkbox"/>
	LIGHTING: fluorescent/LED, general <input type="checkbox"/>	LIGHTING: indirect <input type="checkbox"/>
Security / Nurse Call	NURSE CALL SYSTEM <input checked="" type="checkbox"/>	

Room Fabrics

Room Name | Uptake Room
Room Number UTRM

AusHFG code	Description	Comment
CLTI-002	CEILING: acoustic, drop-in tiles, 600 x 1200	
CLCN-009	CORNICE: wall trim fixing	
DWSC-005	DOOR: 910mm c/o, single leaf, solid core, lead shielding, paint	automated, hold open
DWPR-006	DOOR PROTECTION: protection plate, to 900 AFFL, PVC, prefinished	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-020	WALL FINISH: vinyl, to 1500 AFFL	

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	1	
FQCW-044	BIN: waste, sharps, wall mounted	3	1	
FIHR-004	BRACKET: for display screen, ceiling mounted	1	1	
FIHR-077	BRACKET: for observation camera, ceiling mounted	2	1	
FIHR-016	BRACKET: for sharps bin	2	1	
FIHR-019	BRACKET: for suction bottle	2	1	
ITCL-012	CAMERA: observation	2	1	
MMBE-091	CHAIR: recliner, clinical, electric	3	1	consider bariatric capacity as required
FIDI-001	DISPENSER: alcohol based hand rub	2	1	
FIDI-014	DISPENSER: disposable glove	2	1	
ITBU-012	DISPLAY SCREEN: patient entertainment system	1	1	
MGFP-004	FLOWMETER: oxygen	3	1	
FIHR-029	HOOK: equipment	1	2	
LIFX-027	LIGHT: examination, ceiling mounted	1	1	optional
MMMI-219	PASS THROUGH: injection port	1	1	optional
MMGE-100	SUCTION BOTTLE	3	1	

Engineering Services

AusHFG code	Description	Group	Qty	Comment
ELPR-001	BODY PROTECTED ELECTRICAL AREA	1	1	
ITCL-004	BUTTON: nurse call, emergency	1	1	
ITCL-008	BUTTON: nurse call, staff assist	1	1	
ELGP-208	GPO: double, wall mounted	1	1	
ELGP-122	GPO: emergency power, single, on services panel	1	2	
ELGP-101	GPO: single, ceiling mounted	1	3	
ELGP-109	GPO: single, wall mounted	1	1	
ITCL-028	HANDSET: nurse call, with patient to staff call button	1	1	
ITSE-029	INTERCOM: audio, wall mount	1	1	handsfree to control room
ITIN-011	OUTLET: data, single RJ45, ceiling mounted	1	1	
ITIN-070	OUTLET: MATV	1	1	to patient entertainment screen, may be replaced with data outlet in accordance with ICT policy

Room Name | Uptake Room
Room Number UTRM

AusHFG code	Description	Group	Qty	Comment
MGAS-042	OUTLET: oxygen (O2), on services panel	1	1	
MGAS-062	OUTLET: suction, on services panel	1	1	
MMSP-026	SERVICES PANEL: wall mounted	1	1	
ELBO-012	SWITCH: light, dimmer	1	2	inside and outside room, to be reviewed to suit operational requirements